
Brief on ANSEL Data Acquisition with the DDC-8 DSP

Original by Eryk Druszkiewicz
Edited by Dev A. Khaitan,
Ben Hmiel & WUS

Precaution

ALWAYS inspect analog and trigger signals on the oscilloscope, check relative timing --

BEFORE connecting the cables to the DDC-8 channels !

Signals $\gg 2V$ have the potential of damaging the circuitry of the instrument (a very expensive mistake!)

DDC-8 Software



“BlackBox” programming environment


**Graphical user interface for controlling
the DDC-8DSP**

Preinstalled on all ANSEL computers

Using the BlackBox DDC-8 Interface




1. Start BlackBox ANSEL



2. Open main control panel

Using the BlackBox interface

3. Initialize the connection with the DDC-8DSP board (press only once)



4. You should see a "zero".

If you see a -22 then:

- Close BlackBox
- Reset the USB with the button on the front panel
- using Cypress EZ-USB reprogram the USB chip
- Start-up BlackBox again


If you see a -5 then probably you pressed The Initialize USB button more then once

- close and reopen BlackBox

Using the BlackBox DDC-8 Interface

7. Initialize the Waveform Preview (#Events/Buffer=1)

Waveform Preview



Set the offset for each input channel :

$32 = -1V$ Use 32 for unipolar signals >0

$128 = 0V$

$225 = +1V$

Preview of the Peak values for Ch0 and Ch1

Using the BlackBox DDC-8 Interface

Histogram Digitized Spectrum (Bar Diagram of Frequency vs. Amplitude)

8. Configure desired parameters:


Normally

of Bins: 4096 (4kB)

Offset 2048

Bin(width): 16


9. Initiate histogram



Using the BlackBox DDC-8 Interface

7. Initialize the Waveform preview

Set # of events per Buffer to 1



Waveform values range between:
-32768 to 32764

Control the offset
(30 = -1V, 128=0V, 225=1V)

Input signal + offset => ADC


Preview of the Peak values for Ch0 and Ch1

Preview of the Sum values for Ch0 and Ch1


BlackBox DDC-8 Interface: Experiment NaI(Na-22)

8. Configure desired the parameters

9. Initiate histogram



Using the BlackBox interface



13. Start acquiring data

10. Set number of events to capture, 0 = infinite

11. Configure file outputs

12. Number of samples to define signal and peak. Should be long enough to avoid pileup.

DDC-8/Blackbox Quickstart User's Guide

- 1) Turn on the DDC-8. Wait to hear **3 dings** come from the computer speakers
- 2) Open up blackbox and choose **DDCT>Main** from the top menu
- 3) Begin the initialization procedure
 - a) click on Initialize USB (only once!). You should see "**Interface: 0**" in the log window
If you see something else (possibly -5 or -22), close Blackbox, push the button labeled "Reset USB" on the front panel of the DDC-8 (or power-cycle the device)
 - b) Open the waveform preview window, click "**Initialize Plots**"
Several flat lines should appear, as well as a legend
Double click on the plot to auto-arrange the legend
 - c) Choose the offsets for each channel with the boxes on the right hand side.
For most cases you will want to input "**32**". This will make the offset approx -0.8V, enabling full scale digitization of ~1.8V pulses.
Other offsets can be used here (namely if negative [use "**225**"] or Bipolar [use "**128**"] pulses are to be expected.) For now, stick to an offset of "**32**".
- 4) Inspect the waveform and ensure the DDC-8 is reading the proper signal.
 - a) Attach the signal to be measured to **Ch_0** and the trigger signal to **NIM_IN_0** BNC panels
Remember to inspect them in the scope first to ensure they are <2V!
 - b) Check the box for "**capture waveforms**"
 - c) Choose a value for **# of events/buffer**.
This adjust how often the GUI will update information from the FPGA. **25-50** is sufficient for high count rates. If the event rate is low, you can bring this down to single digits.

DDC-8/Blackbox Quickstart User's Guide

- d) Set each desired Ch#s to **max** and **ena**
- e) Press Apply to save all of the settings adjusted so far.
- f) Adjust the post trigger delay
 - This lets the DDC-8 know after how many samples (40ns ea.) after the trigger it should stop capturing the waveform. **80** is a good value to start with, but you can explore the effect by changing it!
- g) Click "send delay to board".
- h) Click run to see the waveform of the signal. (don't bother saving it)
- 5) After observing that the waveform is correctly displaying, you want to see a histogram of pulse heights
 - a) From the top menu, select **DDCT>Histograms** to display the histogram window
 - b) Configure the histogram parameters, then initialize the histogram
 - Good settings to use are **4096** Bins, **2048** Offset & **16** Bin. These should be the default
 - c) This is just for the initial preview, you will be able to re-bin your data in **Igor** at home
 - d) Back in the main menu, Uncheck "**capture waveforms**" and check "**capture peaks**"
 - ensure that you have set the "**max**" and the "**ena**" of the proper channels, if you indeed want to collect the peak height. The device is too slow to capture waveforms and peak heights simultaneously.
 - e) Click "**Apply**", then "**Send Delay to board**" to collect data and **view the histogram**.
- 6) After ensuring the waveform & histogram are properly displayed, it is time to collect data by checking the "**Write Peaks to File**" box. You may set a specific # of events to collect for a simple energy spectrum (**60,000** - **120,000** should be sufficient depending on the activity of your source and detector used)
 - Or you can put "**0**" to stop collection manually.

Tips/Troubleshooting

- 1) After changing any setting, you need to press "**Apply**", then "**Send delay to board**" before starting data collection. This finalizes the changes in the software, then sends them to the FPGA.
- 2) Be sure to write down in your logbook the Runtime & # of events collected after each run.
 - a. Don't worry if you forget, however, as this information is accessible from the output file.
- 3) Be sure to save all of your data on the EXPERIMENTAL_DATA partitions on the hard drives. Most of the C: drives will fill up quickly.
 - a. This was a common problem in previous years for students. If you data is not saving correctly, be sure to check and make sure the HDD isn't full. your
- 4) Keep an eye on the Log Window. If the machine starts to display "-5"s, there is likely a problem with the connection and its better to power cycle the instrument and reinitialize.
- 5) Sometimes, you will hear the "ding" of the USB disconnect. This happens occasionally via static discharge or jostling of cables. Power cycle the box & reinitialize to continue
- 6) If you want to quickly export the histogram from Blackbox, double click on the legend entry for "**Hist_0**" and a window will appear with the bin values in a list. This makes for easy importation into IGOR, but you will want to get used to doing the histograms manually on your own.
- 7) If you want to adjust the scaling on the histogram to zoom to a particular region, double click on either axis to bring up the menu.
 - a. Be sure to uncheck the "autoscale" box in order to change the scale of a particular axis.

Experimental Setup

- Main Question: What do we want to measure?
 - Is it energy? Is it a timing function? etc.
- Questions to ask before you start:
 - How many events do you need?
 - How long will it take?
 - What could go wrong? How do you know your experiment is proceeding correctly?
- Record keeping is critical for good results. What else should you record?
- All entries should be dated. **NO** entries on loose sheets.
- **ALL** relevant parameters of your experiment and analysis should be recorded in your log book. Do not leave empty pages to complete later. Complete it **NOW**.

Experimental Setup

- What do you record:
 - Type and size of detector and operating condition (e.g. HV setting).
 - Type, location, and strength of sources used.
 - Electronics setup (record model numbers of modules being used).
 - Settings of modules (gain – course and fine, switches).
 - Relative timing of signals, signal shapes (capture scope-screen images).
 - Run conditions (e.g. run length, file name for data, recording mode).
- While you are collecting data, does it make sense?
 - What are checks you can do to ensure data collection is going smoothly?

Experimental Setup

- Signal Chain:
 - PMT -> (50-Ohm term.) Amplifier:
 - Unipolar out -> DDC8 Channel 0
 - Bipolar out -> TSCA -> Gate Generator -> DDC8 Trigger In 0
 - *****At every step** look at the signals. *******
 - Set up the DDC8. Follow **every step** of handy-dandy instruction sheet.
 - Don't cut corners.